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**How to fabricate a CA[®] CLEAR ALIGNER splint using
VECTOR[®] 30 screws and the CA[®] CHROMIUM arch**

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How to fabricate a CA[®] CLEAR ALIGNER splint using VECTOR[®] 30 screws and the CA[®] Arch CHROMIUM

Correction of malocclusion with a CA[®] CLEAR ALIGNER splint gets even more economic, effective and time-saving in comparison to the conventional CA[®] Set Up technique when using the specific VECTOR[®] 30 screw and the CA[®] Arch CHROMIUM.

Relevant literature:
CA[®] Book Volume 2 Advanced
(Chapter: VECTOR[®] 30 Screw with CA[®] Arch CHROMIUM
[English: REF 5309 / German: REF 5307 / Spanish: REF 5306])

Working Material:

REF 2484	VECTOR [®] 30
REF 5317	CA [®] CHROMIUM Arch
REF 3207	ISOFOLAN [®]
REF 3415	DURAN [®] 1.0 mm
REF 3259	DURASPLINT [®] monomer
REF 3258	DURASPLINT [®] polymer
REF 5302	CA [®] Cyano Veneer Fast
REF 3381	DIMO [®] PRO
REF 3226	BIOPERM-Trimmer


Additionally recommended material:

REF 3442	SIL-KITT red
REF 8450	Ortho arch bender
REF 5300	CA [®] pipettes
REF 3220	Blocking out putty
REF 2041	CARAT [®] flat pliers
REF 5372	Set-Up manual saw
REF 5263	CA [®] grinding set
REF 3712	SD polishing brush
REF 5288	Diamond separating disc
REF 3714	SD high gloss polishing disc
REF 3710	SD polishing paste A
REF 3606	SD separating disc, serrated
REF 3460	SD foil scissors A
REF 3461	SD foil scissors B

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Cut off ISOFOLAN® at model rim and insert interdental cuts using a scalpel.

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Embed model into the pellets and thermoform DURAN® 1.0 mm. Then roughly cut it out using the SD foil scissors B.

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
Push the two VECTOR® 30 screws onto the CA® arch CHROMIUM.

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For easier insertion the ends of the CA® arch CHROMIUM are conically tapered.

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
Reduce CA® arch CHROMIUM and bend the retention loops, making sure the retention ends can be used as selective supports on the splint in the molar region.

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If necessary, the lower part of the plastic space holder can be reduced, in order to make sure the screws are positioned in an right angle on the anterior teeth.

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Fix the CA® Arch CHROMIUM and the VECTOR® 30 screws with CA® Cyano Veneer Fast on the DURAN® material. Apply just one drop each to the opening of the bent retentions using the CA® pipette (see fig.)

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Cover selected areas on the splint using blocking-out putty for protection.

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Mix and apply DURASPLINT® according to the instructions and polymerize it in a pressure pot.

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
Then roughly cut out the resin using the CA® carbide bur and mark the saw cuts with a pencil.

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Separate the material within the marked area using the set-up manual saw or separating disc. If necessary, leave the splint on the model during this process for stability reasons.

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Lift off the splint from the model and combine/extend the saw cuts between the screws with the HM tricutter fine (green).

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Shorten the splint using SD foil scissors A and finish it with BIOPERM trimmer.

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A dressed DIMO[®] PRO (see fig.) and SD high-gloss polishing mop in combination with SD polishing paste A are best suited for finishing the rims. Ensure to use them at a low speed (5.000 to 10.000 rpm).

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
Finally conventionally polish the areas where you applied the resin.

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Finished CA[®] CLEAR ALIGNER splint for the oblique segment expansion according to Dr. Pablo Echarri.

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 rotation: 0.35 mm

 rotation: 0.0875 mm

Note: Serious incidents in connection with the use of the product shall be reported immediately to the manufacturer and the competent authority.

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Training courses and seminars.

With the support of our internal and external speakers we offer various courses for advanced training. The training laboratory of our Academy is equipped with cutting-edge technology and can accommodate eight participants. Working in small groups guarantees intense learning combined with individual advice and support.

Our current seminars and training courses for practices and laboratories can be found at:
www.scheu-academy.com

CE

SCHEU-DENTAL GmbH
www.scheu-dental.com

Am Burgberg 20
58642 Iserlohn · Germany

phone +49 2374 9288-0
fax +49 2374 9288-90

SCHEU GROUP

DE/GB 200/01/23 G REF PM0160.01