Occlusion trainer

Due to wearing the protrusion splint, your lower jaw is moved forward during the night. As a result, slight pain in jaw muscles or joints may occur. A temporary bite change is also possible. The occlusion trainer assists you in bringing your jaw back into the normal bite position. To do this, proceed as follows:

- Place one of the enclosed occlusion trainers between your front teeth ensuring that the back teeth are not touching each other.
- While holding the occlusion trainer with your anterior teeth, try to push the lower jaw first forward and then backward.
- 3. Once you moved the jaw to a backward position, try to squeeze the back teeth together.
- 4. Please repeat this procedure until you feel that your teeth are back in the correct bite position.



Hint:

We recommend doing this exercise every morning after removing your splint, to prevent possible bite displacement.

Attention: Serious incidents in connection with the use of the product must be reported to the manufacturer and the competent authority immediately!



Am Burgberg 20 58642 Iserlohn · Germany phone +49 2374 9288-0 fax +49 2374 9288-90 Facebook scheu.dental



