IMPRIMO® LC IBT

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IMPRIMO® LC IBT

Processing Information

| | |
|---------------------------------|---|
| Description: | IMPRIMO® LC IBT is a light-curing resin based on methacrylate and specifically engineered for the fabrication of a transfer matrix in the indirect bonding technique; photo-polymerization is done by means of image projection systems ≤ 405 nm. IMPRIMO® LC IBT has been tested in terms of bio-compatibility and meets both the applicative and mechanical requirements. IMPRIMO® LC IBT should be used in building processes with a z-resolution of 50-100 μm . Printer parameters available for Asiga MAX™ and Asiga PRO 4K. |
| Processing instructions: | Make sure IMPRIMO® LC IBT is tempered prior to use to 23 °C - 30 °C. Fill the material into the tray wearing your personal protective equipment (safety glasses, gloves and face mask). Small air bubbles can be removed carefully with a spatula. In case of contamination of the material such as dirt or fragments, the material can be filtrated and repeatedly mixed due to its low viscosity. Leave the material to stand for approx. one hour to avoid any air inclusions. |
| Contra-indications: | IMPRIMO® LC IBT should not be used for any purpose other than the generative CAD / CAM fabrication of transfer matrixes in the indirect bonding technique. IMPRIMO® LC IBT is a medical product class I. |
| Post-processing: | Upon completion of the printing process, the objects have to be removed from the model platform. Unpolymerized resin residues can be removed completely using isopropanol and a soft brush or may be resolved in an appropriate cleaning device such as IMPRIMO® Clean using IMPRIMO® Cleaning Liquid. |
| Post curing: | Complete curing of the objects is carried out using the light ovens IMPRIMO® Cure or RS cure. The wavelength of 385 nm and the supporting use of protective gas or a vacuum to avoid a dispersion layer is mandatory. The polymerisation in the IMPRIMO® Cure takes 7 minutes with a 20 second nitrogen purge. In the RS cure, polymerisation under vacuum takes 30 minutes. The parameters are already stored in the IMPRIMO® Cure as well as in the RS cure. Deviations from the listed manufacturing process can lead to altered mechanical properties or colour deviations of the material. |
| Shelf life / storage: | Store the material at room temperature (18 °C to 28 °C) and protect it from light. Any exposure to light might damage the material. Carefully seal the container immediately after use. |
| Hazard warnings: | Prolonged skin contact with un-cured material and inhalation of monomer vapours can lead to irritations. In single cases, allergic reactions to certain ingredients of IMPRIMO® LC IBT might occur, for more details, please refer to our safety data sheets on IMPRIMO® LC IBT. In the event of eye contact, rinse thoroughly with water and seek medical advice. In the event of skin contact, rinse thoroughly with running water and soap. Attention: Polymerized resins are chemically resistant, avoid stains on clothes. |
| Disposal: | According to local authority prescriptions. Cured objects can be disposed of with other household waste. Unpolymerized material has to be disposed of at a waste collection point. |
| Hint: | To be used only for the above indication and by qualified personnel. For single use. Serious incidents in connection with the use of the product must be reported immediately to the manufacturer and the competent authority! |
| Technical data: | Density: 1.1 g/ml Viscosity: approx. 2.5 Pa s Cured material: (depends on irradiation unit) Elongation at break: 50 % Tear strength: 6,2 MPa Shore hardness: 40 D Comply: Biocompatibility (irritation and delayed-type hypersensitivity; cytotoxicity) <small>*These data result from measurements of a representative sample.</small> Contains: urethane methacrylates, initiators, stabilizers and dyes. Material Safety Data Sheet: www.scheu-dental.com/en/service |

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