1. Intended use/Indication

Light-curable polymerizable resin intended to be used in conjunction with extra-oral curing light equipment for the fabrication, by additive manufacturing, of dental appliances. **IMPRIMO® LC Temp** is indicated for the fabrication, by additive manufacturing, of temporary anterior dental restorations. **IMPRIMO® LC Temp It** is indicated for the fabrication, by additive manufacturing, of temporary dental restorations, and for the fabrication of preformed denture teeth to be used in a denture.

2. Contraindication

IMPRIMO® LC Temp and **IMPRIMO® LC Temp It** restorations are contraindicated...

- ... if a patient is known to be allergic to any of the ingredients
- ... for bridges with more than one pontic.
- ... for every application that is not part of the indication (see above)

3. Patient Target Group

Persons being treated in the context of a dental procedure.

4. Intended Users

Dentist, dental technician

5. Requirements

Software - information available from:

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Hardware (3D Printing) - information available from:

ASIGA · 2/19-21 Bourke Road, Alexandria, NSW 2015, Australia

Hardware (Post Curing) - information available from:

pro3dure medical (CD-1, CD-2) · Am Burgberg 13, D-58642 Iserlohn, Germany

Additional specifications on www.pro3dure.com

6. Material

IMPRIMO® LC Temp and **IMPRIMO® LC Temp It** consists of functional methacrylic resins and inorganic fillers with particle sizes from 0.4 to 3 microns.

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occlusal 1.5 mm (central fissure)
circumferential 1 mm

8. Material parameters

Radiation penetration depth controllable by exposure time

≥ 100 µm (4 mils)

≥ 125 µm (5 mils)

9. Manufacturing process (fig. 1-10)

- Prepare data (CAD & build preparation).
- Choose process parameters (build style, etc.).
- Transfer prepared data to 3D printing device.
- Prepare 3D printing – shake bottle.
- Fill resin tank of 3D printing device.
- Build the parts.
- Clean parts (with IPA ≥ 97 % or equivalent cleaning agent) approx. 4 min. in an ultrasonic bath or equivalent device – precleaning recommended.
- Dry parts (with compressed air until no residues of IPA or equivalent cleaning agent).
- Post curing (10 min.): inert atmosphere recommended (use adequate light curing device).
- Finish parts.

10. Individualization and cementation

The temporary restorations can be individualized by light curing composite materials. The reduction of the restoration should not exceed 0.3 mm (incisal and vestibular area). Please follow the instructions of the light curing composite manufacturers. For cementation common provisional dental cements (free of Eugenol) and adhesives are suitable. The inner areas of the crown should be conditioned by sandblasting or other commonly used techniques. Please follow the instructions and indications of the corresponding manufacturers.


11. Notice

Please follow the instructions provided by the software manufacturer in terms of parameter settings and design recommendations. Please follow the instructions provided by the hardware manufacturer in terms of parameter settings/printing & post curing recommendations. To avoid detrimental effects on material quality do not expose the liquid material to irradiation under any circumstances. Deviations from the described manufacturing processes or storing conditions may lead to different mechanical and optical properties of the material. Ensure personal protective gear during processing. According to the EU Medical Devices Regulation, users/patients are obliged to report serious events with a medical device to the manufacturer and to the competent authority of the country in which they occurred. **Caution:** Avoid breathing dust/fume/gas/mist/vapours/spray. Wash hands thoroughly after handling. Use only outdoors or in a well-ventilated area. Contaminated work clothing should not be allowed out of the workplace. Avoid release to the environment. Wear protective gloves/protective clothing/eye protection/face protection. IF ON SKIN: Wash with plenty of soap and water. IF INHALED: Remove person to fresh air and keep comfortable for breathing. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Call a POISON CENTER or doctor/physician if you feel unwell. Specific advice (see information on this label). If skin irritation occurs: Get medical advice/attention. If skin irritation or rash occurs: Get medical advice/attention. If eye irritation persists: Get medical advice/attention. Take off contaminated clothing and wash it before reuse. Collect spillage. Store in a well-ventilated place. Keep container tightly closed. Store locked up. Dispose of contents/container according to official regulations.

12. Hazard statements

Causes skin irritation. May cause an allergic skin reaction. Causes serious eye irritation. May cause respiratory irritation. Toxic to aquatic life with long lasting effects.

Instructions for use · Gebrauchsanweisung Mode d'emploi · Instrucciones de uso Istruzioni per l'uso

Physical properties*/
Physikalische Eigenschaften*/
Propriétés physiques*/
Propiedades físicas*/
Proprietà fisiche*:

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ASIGA - /2/19-21 Bourke Road, Alexandria, NSW 2015, Australia

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exocad GmbH - Julius-Reiber-Straße 37, 64293 Darmstadt, Germany