

Splint Remover

Using the Splint Remover, the patient can easily remove thermoformed aligners or splints. Please proceed as follows:



1. Hook the Splint Remover into the rear edge of the aligner or splint. Then begin with loosening the splint at one side with a slight downwards movement.
2. Proceed in the same way at the other side.
3. Repeat these steps (alternating between right and left side) until the aligner/splint can be easily removed in order to make sure not to damage it.

Non-sterile.

Ref: 5359.1 / 5 pcs.



SCHEU-DENTAL GmbH
www.scheu-dental.com

Am Burgberg 20
58642 Iserlohn · Germany

phone +49 2374 9288-0
fax +49 2374 9288-90

 Facebook
scheu.dental

SCHEU

DE/GB 5.000/09/23 G REF PM0175.01

