IMPRIMO® LC MJF

Produktinformation

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IMPRIMO® LC MJF

Product information

Description:	IMPRIMO® LC MJF is a light-curing resin based on methacrylate and specifically engineered for the fabrication of a devices for oral and maxillo-facial surgery as well as for implantology; photo- polymerization is done by means of image projection systems (385 nm). Suitable for current sterilization methods such as plasma sterilization, gamma sterilization, ethylene oxide sterilization and autoclave sterilization. IMPRIMO® LC MJF has been tested in terms of bio-compatibility and meets both the applicative and mechanical requirements. IMPRIMO® LC MJF should be used in building processes with a z-resolution of 25-100 µm. Printer parameters available for Asiga MAX™ and Asiga PRO 4K.
Processing instructions:	Make sure IMPRIMO® LC MJF is tempered prior to use to 23 °C - 30 °C. Fill the material into the tray wearing your personal protective equipment (safety glasses, gloves and face mask). Small air bubbles can be removed carefully with a spatula. In case of contamination of the material such as dirt or fragments, the material can be filtrated and repeatedly mixed due to its low viscosity. Leave the material to stand for approx. one hour to avoid any air inclusions.
Contra-indications:	IMPRIMO® LC MJF should not be used for any purpose other than the generative CAD / CAM fabrication of devices for oral and maxillo-facial surgery as well as implantology. IMPRIMO® LC MJF is a medical product class IIa.
Post-processing:	Upon completion of the printing process, the objects have to be removed from the model platform. Unpolymerized resin residues can be removed completely using isopropanol and a soft brush or may be resolved in the cleaning device such as RS wash or IMPRIMO® Clean using IMPRIMO® Cleaning Liquid.
Post curing:	Final and complete curing of the objects is done in the curing unit such as RS cure or IMPRIMO® Cure. The use of a vacuum or protective gas to prevent the creation of a dispersion layer is crucial. A wavelength of 385 nm is required. The polymerization process in the light oven RS cure in a vacuum takes 6 minutes. The polymerisation process in a nitrogen environment in the IMPRIMO® Cure takes 5 minutes. The corresponding parameters are stored in the internal programme of the RS cure unit. Any deviation from the described post-exposure process can lead to modifications of the mechanical properties or the material colour.
Shelf life / storage:	Store the material at room temperature (18 °C to 28 °C) and protect it from light. Any exposure to light might damage the material. Carefully seal the container immediately after use.
Hazard warnings:	Prolonged skin contact with un-cured material and inhalation of monomer vapours can lead to irritations. In single cases, allergic reactions to certain ingredients of IMPRIMO® LC MJF might occur, for more details, please refer to our safety data sheets on IMPRIMO® LC MJF. In the event of eye contact, rinse thoroughly with water and seek medical advice. In the event of skin contact, rinse thoroughly with running water and soap.
Disposal:	Attention: Polymerized resins are chemically resistant, avoid stains on clothes. According to local authority prescriptions. Cured objects can be disposed of with other household waste. Unpolymerized material has to be disposed of at a waste collection point. All information on processing our material - verbal, written or practical - is given to the best of our knowledge and must be read as references without obligation. Material use and processing is beyond our control and lies within the responsibility of the user.
Technical data:	For professional use only. Density: 1.1 g/ml Viscosity: approx. 0.7 Pa s Cured material (depends on irradiation unit): Flexural strength: 117 MPa Elongation at break: 5 % Solubility: 0.51 µg mm-3 Water absorption: 19.9 µg mm-3 Shore hardness: 85 D Comply: biocompatibility (hemocompatibility); irritation and delayed-type hypersensitivity; genotoxicity, carcinogenicity and reproductive toxicity; systemic toxicity; cytotoxicity). *These data result from measurements of a representative sample. Contains: alkoxyated bisphenol-A dimethacrylate, initiators, stabilizers and dyes. Material Safety Data Sheet: www.scheu-dental.com/en/service

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SCHEU-DENTAL GmbH
SCHAU-DENTAL GmbH