

Product Data Sheet

- Pliers - CA[®] TIP -

CA [®] TIP 1	CA [®] TIP 2	CA [®] TIP 3	CA [®] TIP 4
REF 5281	REF 5291	REF 5292	REF 5295

Material:

AISI 420 (stainless steel)

Alloy composition AISI 420 (stainless steel)

C	≥	0,15 %
Si	≤	1,00 %
Mn	≤	1,00 %
P	≤	0,040 %
S	≤	0,030 %
Cr	=	12,00 - 14,00 %
Mo	=	---
Ni	=	---
Other	=	---
Fe	=	Residue