



Laboratório Dentário Larissa Merth  
Rua Conselheiro Joaquim Machado, 57  
8600-682 Lagos/Portugal  
Tlm. 913 223 044  
Tel. 282 094 611  
aligner.algarve@gmail.com  
www.dentalmerth.com



DLT Dental Labor Trunzer GmbH  
Am Ring 1  
74722 Buchen/Germany  
Tel. +49 (0) 6281 - 2450  
aligner.algarve@gmail.com  
www.dlt-trunzer.de



Practice stamp/Praxisstempel



## With confidence an attractive smile.

Made of clear-transparent material, CA® CLEAR ALIGNER is almost invisible, and you can speak and smile during the treatment period without any difficulty. This guarantees that your mouth will be comfortable - almost like normal.

You'll notice when you wear them: the aligners are pleasantly smooth and comfortable, because they don't require brackets or wires. When you eat and brush your teeth or the aligners, simply take CA® CLEAR ALIGNER out.

## The process for ideal treatment results.

### 1 Consultation

Talk with your dentist about whether the CA® CLEAR ALIGNER system is the right treatment for you.

### 2 Impression

The impressions taken from your upper and lower jaws form the exact basis to create the aligners.

### 3 Creating the aligners

Computer-controlled production of the aligners takes place in the lab, precisely made just for you.

### 4 Treatment

Your dentist will monitor treatment results regularly while you wear the aligners. For best results, additional impressions of your teeth are taken and new aligners are made.

**CA** CLEAR ALIGNER



## Beautifies your smile at any time.

The invisible way of teeth correction.

Learn more about  
Katharina's successful  
CA® CLEAR ALIGNER treatment:



Find out more at:  
[www.ca-clear-aligner.com](http://www.ca-clear-aligner.com)

GB 5.000/01/24 G REF 0159.02





# Gentle and contemporary.

Aligning teeth with aesthetic high-quality splints is an important therapy in modern orthodontics. Especially for treating adults, this concept has proven successful. It provides optimal treatment results along with a very high level of comfort.

There is a wide range of applications for the innovative CA® CLEAR ALIGNER splint:



# Your advantages at a glance.

## effective

Treatment will provide optimal, precise results in a short span of time.

## hygienic

We recommend CETRON® cleansing powder to care for this product. Please ask your pharmacist.



## removable

The aligners can be taken in and out of your mouth quickly and simply at any time (e.g. when you're eating and brushing your teeth).

## cost-effective

In many cases, therapy is more affordable than other treatment options for teeth alignment.



## transparent

Almost invisible, you can speak and laugh without problems.

## comfortable




The aligners are pleasantly smooth, without brackets or wires.

# The 3-Splint-System

Three material thicknesses, worn during 4 weeks for a unique and beautiful smile: With each treatment step, you'll be wearing 3 splints of different thicknesses.

Thanks to this steady increase in pressure, your teeth are aligned efficiently, yet smoothly and with extraordinary comfort.

For easy distinction, the corresponding degree of hardness is marked on the aligners in the posterior segment.

STEP 1	Week 1	soft	
	Week 2	medium	
	Week 3	hard	
	Week 4	hard	
STEP 2	Week 5	soft	
	Week 6	medium	
	Week 7	hard	
	Week 8	hard	
STEP 3	Week 9	soft	
	Week 10	medium	
	Week 11	hard	
	Week 12	hard	



Ask your doctor about treatment and he or she will be happy to advise you.

[www.ca-clear-aligner.com](http://www.ca-clear-aligner.com)